SYNTROPIQ
SPINAL
IMPLANTS

TECNICA QUIRÚRGICA

SURGICAL TECHNIQUE

PERSEUS

PLIF

DYNAM'X
SPINAL INTERBODY CAGE SYSTEM
SISTEMA CAJAS INTERSOMATICAS

PORTA IMPLANTES

The universe is ruled by Entropy,
the unforgiving tendency of decay
into disorder. Domination of entropy
in a life-less universe is inevitable.

It is the life that makes the difference.
Life is the moving force that creates
the complex out of the basic, **Syntropy**.

Syntropy is the tendency towards
energy and matter **concentration**,
tendency to create order out of disorder.

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Syntropy, noun
the tendency towards energy concentration,
order, organization and life.

From Greek **syn** "together" + **tropos** "tendency".
The term was first coined by the mathematician
Luigi Fantappiè in 1941 as a term contrary to entropy.

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ABOUT US

We are a privately held spine interbody implant and surface technology company. Our team members with rich experience more than 15 years in spine have seen over 10 000 spine surgeries and had hundreds of hours of discussions with you – spine surgeons – to obtain much appreciated feedback on how to treat patients with spine disorders most successfully. And here are cool things you are going to like about our implants:

TITANIUM POWDER MATERIAL

- one of the best materials for bone fusion in the medical industry.
- its low density allows MRI and CT diagnostic compatibility.
- highly porous implants with a large contact area to achieve the best fusion.

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* RESPECTING SAGITTAL BALANCE

- to avoid sagittal imbalance in spine our broad product portfolio offers different angulations of the cages including hyperlordotic cages.

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

* DYNAMIC CONCEPT

- our Dynam'X Interbody Cage System respects Wolff's Law.
- the cage's design allows dynamization of fusion by controlled elastic microdeformation of the contact net. This phenomenon stimulates dynamization of the bone rebuilding process in accordance with Wolff's Law.
- the open construction of the implant facilitates free blood flow and acceleration of bone formation to achieve best fusion.

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

Items highlighted in grey are available upon request.

IMPLANTS

IMPLANTES

10×21 (w × l) Anatomic, Hollowed, XSmall

h	5°	8°	12°
7 mm	2410210507	2410210807	2410211207
8 mm	2410210508	2410210808	2410211208
9 mm	2410210509	2410210809	2410211209
10 mm	2410210510	2410210810	2410211210
11 mm	2410210511	2410210811	2410211211
12 mm	2410210512	2410210812	2410211212
13 mm	2410210513	2410210813	2410211213
14 mm	2410210514	2410210814	2410211214

10×25 (w × l) Anatomic, Hollowed, Small

h	0°	5°	8°
7 mm	2410250007	2410250507	2410250807
8 mm	2410250008	2410250508	2410250808
9 mm	2410250009	2410250509	2410250809
10 mm	2410250010	2410250510	2410250810
11 mm	2410250011	2410250511	2410250811
12 mm	2410250012	2410250512	2410250812
13 mm	2410250013	2410250513	2410250813
14 mm	2410250014	2410250514	2410250814

10×29 (w × l) Anatomic, Hollowed, Medium

h	0°	5°	8°
7 mm	2410290007	2410290507	2410290807
8 mm	2410290008	2410290508	2410290808
9 mm	2410290009	2410290509	2410290809
10 mm	2410290010	2410290510	2410290810
11 mm	2410290011	2410290511	2410290811
12 mm	2410290012	2410290512	2410290812
13 mm	2410290013	2410290513	2410290813
14 mm	2410290014	2410290514	2410290814

10×33 (w × l) Anatomic, Hollowed, Large

h	0°
7 mm	2410330007
8 mm	2410330008
9 mm	2410330009
10 mm	2410330010
11 mm	2410330011
12 mm	2410330012
13 mm	2410330013
14 mm	2410330014

Items highlighted in grey are available upon request.

INSTRUMENTS

INSTRUMENTAL

IMPLANT TRIAL

height	25 mm, 5°	33 mm, 0°
7 mm	202507	203507
8 mm	202508	203508
9 mm	202509	203509
10 mm	202510	203510
11 mm	202511	203511
12 mm	202512	203512
13 mm	202513	203513
14 mm	202514	203514

INTERBODY DISTRACTOR

height	25 mm	33 mm
7 mm	201007	211107
8 mm	201008	211108
9 mm	201009	211109
10 mm	201010	211110
11 mm	201011	211111
12 mm	201012	211112
13 mm	201013	211113
14 mm	201014	211114

DISC SHAVER

height	25 mm	33 mm
7 mm	220007	230007
8 mm	220008	230008
9 mm	220009	230009
10 mm	220010	230010
11 mm	220011	230011
12 mm	220012	230012
13 mm	220013	230013
14 mm	220014	230014

DISC SHAVER – QUATRO

height	25 mm	33 mm
7 mm	221007	231007
8 mm	221008	231008
9 mm	221009	231009
10 mm	221010	231010
11 mm	221011	231011
12 mm	221012	231012
13 mm	221013	231013
14 mm	221014	231014

INSTRUMENTS

INSTRUMENTAL

DESCRIPTION	REF NUMBER	
IMPLANT HOLDER PORTA IMPLANTES	200001	
SHAFT FOR IMPLANT HOLDER MANGO CON ROSCA LIBRE IMPACTO	200002	
IMPLANT HOLDER – SIDE OFFSET PORTA IMPLANTES LATERAL	200003	
IMPLANT HOLDER WITH HAMMERING ZONE PORTA IMPLANTES CON ZONA DE MARTILLO	200004	
HANDLE – T-SHAPE MANGO EN T	200005	
HANDLE – STRAIGHT MANGO RECTO	200006	
HANDLE – STRAIGHT FOR HAMMERING MANGO RECTO PARA MARTILLAR	200006A	
ENDPLATE CURETTE – RECTANGULAR RASPADOR PLATILLO RECTANGULAR	200055	
ENDPLATE CURETTE RASPADOR PLATILLO	200056	
ENTRANCE CHISEL – 6×8 MM ESCOPLO 6 X 8 MM	200057	
ENTRANCE CHISEL – 8×10 MM ESCOPLO 8 X 10 MM	200058	
HAMMER MARTILLO	200059	
IMPLANT FILLING SUPPORT / Soporte Caja	200060	
TRIAL EXPACTOR EXTRACTOR DE PRUEBA	200061	
IMPLANT REMOVAL TOOL* EXTRACCION DE IMPLANTES	199100	

*) Not in any standard set configuration.
Delivered separately upon special
request.

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
PATIENT POSITIONING

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
The position of the patient shall expose the spine level to be fused. Surgeons shall evaluate the most appropriate position taking into consideration approach, decompression procedure and fusion technique. By using C-Arm proper level of lumbar spine must be verified. The incision should allow the adequate approach to reach the targeted spine segment(s). Additional instruments like vertebral distractor, soft tissue retractors are subject to additional room requirements.

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In order to have a comfortable surgical field, tissue retractor system is highly recommended. It is up to surgeon experience to define and perform soft tissue approach and bone decompression procedure. In some cases specific patients positioning may be required.

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
DISC REMOVAL AND ENDPLATE PREPARATION

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1. Standard laminectomy for classical PLIF approach should be performed

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2. Syntropiq Chisel allows to create open window for discectomy procedure. It comes in two sizes 6x8 mm and 8x10mm.

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
3. Syntropiq offers different types of endplate scrapers to facilitate the preparation of the endplates. The appropriate choice depends on the approach, surgical preferences or expected results.

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4. The discectomy procedure should be performed using standard instruments. **Perseus Lumbar Cage** incorporates titanium net which allows free flow of bone cells. Rasp and scrapers are recommended when preparing the endplate in order to create efficient bone contact.

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
IMPLANT TRIALING

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1. In order to select the appropriate size of the **Perseus** cage implant, several trials of different sizes of implants shall be used. Attach implant trial to the straight handle with hammering zone. By slightly hammering insert the implant trial into intervertebral disc space, increase the high of trial until firmly seat in intervertebral disc space will be achieved. Keep the same trajectory and angle as intended for the introduction of the implant.

2. While the trial is performed, lateral x-ray is highly recommended in order to assess the high, angulations and length of the implant. In case of difficulties of removing the trial Syntropiq expactor may be used.

Please note the trials are available in 25 mm length, 5 degrees and 33 mm, 0 degrees. Please note that implants are available in 0, 5, 8 and 12 degrees of lordosis. If larger lordosis restoration is intended, please select cage with 8 or 12 degrees gradient.

3. Syntropiq offers 3 types of handles:
T-Shape – dedicated to rotational tools, like scrapers
Straight – intended for insertion, like trials
Straight with hammering – intended for insertion when hammering forces are required

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IMPLANT PREPARATION

PREPARACION DE LA IMPLANTACION


1. Perseus Cage is intended to be filled by bone substitute or harvest allograft. Both substitution can be delivered into the cage. Implant can be filled before or after implantation procedure. In case of postimplantation filling procedure only paste or gel formulation can be used.

2. Please place the cage into the Filling Support.
By using Syringe inject the bone substitute into the cage.

PORT

3. Perseus System offers unique possibility to fill cage by autograft. Model Hollowed has incorporated open net which facilitates grafting. Place cage with hole up and use standard instruments to place and push bone graft into the cage. For hollowed version availability please contact with your local Syntropiq Representative.

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4. Perseus System offers two implant inserters.

Regular Implant holder with comfortable handle for ordinary application and side offset for limited eyes view or microscope or MIS application. Connect the selected Implant Holder to Cage and secure the position by clockwise turning the knob on silicone handle.

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IMPLANTATION

1. Increase distraction to facilitate the introduction of the device. Introduce the cage hammering gently. Keep appropriate angle to introduce the device. Lateral x-ray is highly recommended during the implantation procedure. After confirming the final position, detach the Implant Holder by turning the knob anticlockwise. Remove the Implant Holder and verify the implant position both in AP and lateral x-rays.

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2. Before placing the second implant, autologous bone or bone substitutes should be injected into the intervertebral disc space.

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Please repeat the procedure described above for the second cage on the other side.

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
3. The **Perseus Cage** can be filled after implantation using cannula system. Please ask Syntropiq Representative for more information.

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IMPLANT REMOVAL

1. Attach either the Implant Holder (Inserter) or universal Implants Removal Tool in a clockwise rotation to the implant. Be careful to avoid pushing the implant anteriorly. A slap hammer or slotted mallet may be used in conjunction with the Inserter or Implant Removal Tool for removal of the implant if desired.

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To use, apply an upward force to the slap hammer. Repeat until implant is removed from the intervertebral disc space.

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2. If distraction was utilized during implantation, be sure to re-apply distraction to allow easier removal of the implant. Vertebral bone overgrowth or osteophytes may be removed to facilitate retrieval of the implant.

Notes:

An osteotome can be used at the interface between the Implant and endplates to disengage the construct. Use of distraction is suggested to allow easier access to the implant/endplate interface.

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