


Distribuidor por:



BIOSpine. S.L.

Camino de los Cinamomos 283-5
33203 - Gijón (Asturias)
985195325
info@biospine.info

MEDICAL
adi+us[®]



Cemad

**Vertebral Radiopaque
Bone Cement**

Cemad Vertebral Radiopaque Bone Cement

Cemad is a radiopaque, low viscosity and acrylic based (polymethylmethacrylate) bone cement. It contains a pouch with radiopaque powder sterilized by ethylene oxide and an ampoule with liquid sterilized by filtration.

Key Features

Reduced mixing time

Unique dynamic mechanical properties


Excellent handling characteristics

Low viscosity

Ideal for minimally invasive surgery


ISO 5833 standard

Working Time According to ISO 5833



Cemad Bone Cements are CE marked and complying ISO 5833 specifications.

Mechanical Properties According to ISO 5833



Order Information



Reference Code	Product Description	Contents
ADCEMLV	CEMAD-LV Vertebral Radiopaque Bone Cement	20g
ADCEMLVF	CEMAD-LVF Vertebral Radiopaque Bone Cement	20g