

ALAYA SPINE FIDUCIAL – DATA SHEET

General Description

The ALAYA SPINE Fiducial is intended for use with ALAYA SPINE Robotic Assistant as an aid during image registration process to perform the recording of the patient's position by intraoperative imaging.

It provides an accurate robotic surgical solution to people with spinal pathologies (vertebral fractures, spinal stenosis, and kyphosis, among others) who need a surgical intervention of transpedicular screw insertion.

Regulatory Compliance

The product meets safety and performance requirements of MDR (EU) 2017/745 for Class I sterile medical devices.

CYBER SURGERY S.L.

Paseo Mikeletegi 71, 3º Local 336
Donostia-San Sebastián 20009, Gipuzkoa (Spain)
Tel : +34 943 12 46 27
Email : customerservice@cyber-surgery.com

Technical Details

Material: Polyether ether ketone (PEEK)

Spheres material: Stainless Steel 316L

Kit size: 200 x 180 x 40 mm

Sterilization method: Gamma radiation

Primary packaging: Tyvek pouch

Secondary packaging: Cardboard box

Weight: 0.03 kg

Shipping and storage temperature and humidity:
Under 25°C / 30 to 65%, without condensation.

